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
Table 1. Colour reaction of Narcotic Drugs


Reagent	Marquis	Fröhde	Mecke
Sample	Reddishviolet	Reddishviolet ~Darkgreen	Green~Blue ~Brownishgreen
Heroin	Deepviolet	Violet~ Greenishbrown	Deepblue~ Greenishbrown
Acetylcodein	Violet	Darkgreen~ Brownishgreen	Deepviolet ~ ~Bluegreen ~Brown
Cocaine	Colourless	Colourless	Colourless

(1) (2) (3) (4)

Fig.1.Thin-layer Chromatograms of Narcotic Drugs

- (1) Sample
(2) Heroin
(3) Acetylcodeine
(4) Cocaine
Solvent : Chloroform-Acetone Methanol(5 : 4 : 1)

(1) (2) (3)

Fig.2.Thin-layer Chromatograms of Narcotic Drugs.

- (1) Hydrolysis of Sample
(2) Morphine
(3) Codeine
Solvent : Chloroform-Methanol(4 : 1)

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Table 2. Colour reaction of Narcotic Drugs

Reagent	Marquis	Fröhde	Mecke
Sample	Reddishviolet	Brownishviolet ~Greenishbrown	Greenishbrown ~Green
Morphine	— 〃 —	Violet~ Greenishbrown	Deepblue~ Greenishbrown
Ethyl-Morphine	Deepviolet	Yellowishgreen ~Blue	Deepblue~ Green~Brown


Fig. 3. Thin-layer Chromatograms of Narcotic Drugs.


- (1)Sample
(2)Morphine
(3)Ethylmorphine
(4)Codeine
Solvent : Chloroform-Methanol
(4 : 1)

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Table 3. Colour reaction of Narcotic Drugs


Reagent	Marquis	Fröhde	Mecke
Sample	Deepviolet	Violet~ Greenishbrown	Blue~Green ~Brown
Heroin	— 〃 —	— 〃 —	— 〃 —

(1) (2)

Fig. 4. Thin-layer Chromatograms of Narcotic Drugs.


(1) Sample
(2) Heroin
Solvent : Chloroform-Aceton-Methanol(5 : 4 : 1)

(1) (2)

Fig. 5. Thin-layer Chromatograms of Narcotic Drugs


(1) Hydrolysis of Sample
(2) Hydrolysis of Heroin
Solvent : Chloroform-Aceton-Methanol(5 : 4 : 1)

(1) (2)

Fig. 6. Thin-layer Chromatograms of Narcotic Drugs

(1) Sample
(2) Powdered Opium
Solvent : Chloroform-Aceton-Methanol(5 : 4 : 1)
a : Morphine
b : Codeine
c : Thebaine
d : Papaverine
e : Noscapine

(1) (2)

Fig. 7. Thin-layer Chromatograms of Narcotic Drugs.

(1) Sample
(2) Powdered opium
Solvent : Aceton-Methanol (4 : 1)

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Table 4. Colour reaction of Narcotic Drugs

Reagent	Marquis	Fröhde	Mecke
Crude Morphine	Reddishviolet	Peddishviolet ~ Green ~ Brownishgreen	Green ~ Greenishbrown
Powdered Opium	—	Violet ~ Brownishgreen	Green ~ Greenishbrown

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Identification of Narcotic Drugs.

Mitsuko YOSHIDA
YOKOHAMA Customs Laboratory.
1-1 Kaigan-dori Nakaku Yokohama city.

Received Sep.29,1969