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Sample concentration : 5 mg/ml

Injection volume : 1  $\mu$ l


Fig. 2 Effect of Water contents on the separation of aromatic sulfonates

Statfonary phase:  $\mu$ -Bondapak NH<sub>2</sub>  
The concentration of PIC - A reagent in each  
Mobile phase is  $2.5 \times 10^{-3}$ M  
pH:7.0  
Flow rate:1.0ml/min

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Table 1 Effect of the type of counter ion on the separation of aromatic sulfonates

counter ion	Elution time (min) of compound			none
	B. M. S	N. M. S ( $\alpha$ -)	N. D. S (2,6-)	
Tetra ethyl ammonium chloride	5.3	7.2	14.8	
Tetra propyl ammonium chloride	5.8	7.6	15.2	
Tetra butyl ammonium chloride	6.2	7.8	15.8	
PIC-A reagent	6.2	8.0	16.2	

Stationary phase:  $\mu$ -Bondapak NH<sub>2</sub> 4mm × 30cm

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
Fig. 3 Effect of molar concentration of counter ion in mobile phase on the separation of aromatic sulfonates

Stationary phase:  $\mu$ -Bondapak NH<sub>2</sub>

Mobile phase: Water: Methanol=(80:20 v/v)

adjusted to pH 7.5 with acetic acid


Fig. 4 Effect of pH of mobile phase on the separation of aromatic sulfonates

1. B.M.S
  2. -naphthalene sulfonate (-N.M.S)
  3. -N.M.S
  4. 2.6-N.D.S
  5. 1.5-naphthalene disulfonate (1.5-N.D.S)
  6. 2.7-naphthalene disulfonate (2.7-N.D.S)
- Stationary phase;  $\mu$ -Bondapak NH<sub>2</sub>  
Mobile phase; Water : Methanol = (80:20 v/v)  
containing  $2.5 \times 10^{-3}$ M  
PIC-A reagent

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Table 2 The recommended conditions for the separation of aromatic sulfonates.

Stationary phase	$\mu$ -Bondapak NH <sub>2</sub> , 4 mm × 30cm
Mobile phase	Water : Methanol = 80 : 20 (v/v) containing $2.5 \times 10^{-3}$ M PIC-A reagent
pH of mobile phase	7.0
Flow rate	1 mL/min
Sample concentra-tion	5 mg/mL
Injection volume	1 $\mu$ L
Detector	UV detector (254nm)
Temperature	Room temperature


Fig. 5 H.P.L chromatogram of aromatic sulfonates

Stationary phase;  $\mu$ -Bondapak NH<sub>2</sub>  
4mmID 30cm  
Mobile phase; Water : Methanol = 80:20 (v/v)  
containing  $2.5 \times 10^{-3}$ M PIC-A  
reagent  
pH of mobile phase; 7.0  
Flow rate; 1.0 mL/min  
Sample concentration; 5 mg/ml  
Injection volume; 1  $\mu$ L

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**Separation of Aromatic Sulfonates by High performance Liquid Chromatography  
Organooionic Compounds by Ion - pair Chromatographic Techniques ( 3 )****Separation of**

Hironobu OSHIRO , Soei SATO and yukio OHNO\*

\*Central Customs Laboratory , Ministry of Finance ,

531 , Iwase , Matsudo - Shi , Chiba - Ken , 271 Japan

The conditions for the separation of benzene monosulfonate , 1,3-disulfonate , 1,5-disulfonate and naphthalene monosulfonates , three isomers of naphthalene disulfonate were investigated.

The retention and separability of aromatic sulfonates were drastically affected with pH , amount of water and molar concentration of counter ion in mobil phase

The aromatic sulfonates were separated by following conditions.

Column : μ - Bondapak NH<sub>2</sub> (10 μ of amino propyl silan)

packed in 4mmID×30cm stainless column

Mobile phase : Water : Methanol ( 80:20v/v )

adjusted to pH 7.0 with acetic acid ,

containing  $2.5 \times 10^{-3}$ M PIC - A reagent

Detector : UV ( 254nm )

Sample concentration : 5mg/ml

Flow rate : 1ml/min

Injection volume : 1 μl ( 5 μg )

- Received Sep . 16 , 1980 -