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Table 1. Samples of Tannin Extract

Tannin Extract	Name of Samples	Origin	Sample №	Appearance of Samples		Remarks (received from)
				Colour	State	
Chesnut	Chesnut Extract Powder	Italy	1	Maroon	Powder	Yokohama Customs Lab.
	Chesnut Extract Dulcochrom	Italy	2	"	"	"
	Chesnut Extract Dulcotan	Italy	3	"	"	"
	Chesnut Extract Rit an	Italy	4	"	"	"
Gambier	Gambier Extract	Singapore	5	Dark Brown	Small Cubic	"
	Gambier Extract Powder	Singapore	6	Brown	Block Powder	"
Heamlock	Heamlock Extract Powder NDK	?	7	Reddish Brown	Powder	"
Mangrove	Mangrove Extract	New Guinea	8	Black	Glistening Amorphous Solid	"
	Mangrove Extract R.Brand	Borneo	9	"	"	"

Tannin Extract	Name of Samples	Origin	Sample №	Appearance of Samples		Remarks (received from)
				Colour	State	
Myrobalan	Myrobalan Extract	India	10	Black	Glistening Amorphous Solid	Yokohama Customs Lab.
	Myrobalan Nut (not Extract)	India	11	Yellowish Brown	"Nut"	"
	Myrobalan Powder	U.K.	12	Yellowish Brown	Powder	"
Quebracho	Quebracho Extract (Soluble) Crown Brand	Argentina	13	Dark Brown	Glistening Amorphous Solid	"
	Quebracho Extract (Insoluble) C. F. del C. Brand	Argentina	14	Brown	"	"
Sumac	Sumac Extract Powder	Italy	15	Yellowish Brown	Powder	"
Wattle (Mimosa)	Wattle (Mimosa) Bark (not Extract)	Africa	16	Light Brown	"Bark"	"
	Wattle (Mimosa) Extract Elephant Brand	South Africa	17	Black	Glistening Amorphous Solid	"
	Wattle (Mimosa) Extract Kenmosa Brand	East Africa	18	"	"	"
	Wattle (Mimosa) Extract Rhomosa Brand	Rhodesia	19	"	"	"
	Wattle (Mimosa) Extract	South Africa	20	"	"	Tokyo Customs Laboratory
?	Hausruk Tannin	Holland	21	Black	Glistening Amorphous Solid	Yokohama Customs Lab.
?	Myrtan	Australia	22	"	"	"
-	Tannic Acid (Japan Pharmacopeia)	-	23	Yellowish Brown	Amorphous Powder	-

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Table 2. Appearance of Extract Solutions

Extraction Solvent	Tannin Extract	Sample №	Colour of Extract Solution	Colour of Extract Solution under U.V. Light
Ethyl-acetate (from sodium bicarbonate solution)	Chesnut	1	Light yellowish brown	Lavender
	Gambier	5	Reddish brown	Bright yellow(^{slightly} greenish)
	Mangrove	8	Light Orange	Rose-pink
	Myrobaran	10	Light brown	Dull yellow
	Quebracho	13	Red	Bright yellow
	Whattle	17	Red	Greenish yellow
50% Methanol	Chesnut	1	Dark brown	Pale yellowish green
	Gambier	5	Reddish brown	Pale yellowish green
	Mangrove	8	Dark reddish brown	Pale pink
	Myrobaran	10	Brown	Yellowish green
	Quebracho	13	Reddish brown	Yellowish green
	Wattle	17	Reddish brown	Pale yellowish green
50% Acetone	Chesnut	1	Dark brown	Pale yellowish green
	Gambier	5	Dark reddish brown	Pale yellowish green
	Mangrove	8	Dark reddish brown	Pale reddish orange
	Myrobaran	10	Dark brown	Yellowish green
	Quebracho	13	Dark reddish brown	Yellowish green
	Wattle	17	Reddish brown	Yellowish green
Ethyl-ether n-Butanol Aceton (2:3:3)	Chesnut	1	Pale yellow	Lavender
	Gambier	5	Yellow	Bright yellow(^{slightly} greenish)
	Mangrove	8	Pale yellow	Rose-pink
	Myrobaran	10	Light yellow	Dull yellow
	Quebracho	13	Light reddish orange	Bright yellow
	Wattle	17	Light yellowish orange	Greenish yellow
Ethyl-acetate (from sodium sulfite solution)	Chesnut	1	Light yellowish brown	Lavender
	Gambier	5	Reddish brown	Bright yellow(^{slightly} greenish)
	Mangrove	8	Light orange	Rose-pink
	Myrobaran	10	Light brown	Dull yellow
	Quebracho	13	Red	Bright yellow
	Wattle	17	Red	Greenish yellow

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Table 3. Appearance of Ethyl - acetate Solutions under Ultraviolet Light

Tannin Extract	Sample №	Colour	Tannin Extract	Sample №	Colour
Chesnut	1	Lavender	Quebracho	13	Bright yellow
	2	Lavender		14	Bright yellow
	3	Lavender	Sumac	15	Dull orange
	4	Lavender		16	Dull yellow
Gambier	5	Bright yellow (slightly greenish)	Wattle	17	Greenish yellow
	6	Bright yellow		18	Light greenish yellow
Heamlock	7	Pale blue		19	Light greenish yellow
	8	Pale blue		20	Light greenish yellow
Mangrove	9	Rose-pink		21	Violet
	10	Rose-pink~violet		22	Lavender
Myrobalan	11	Dull yellow (slightly greenish)	(Hausruk Tan.)	23	Pale violet
	12	Dull yellow (")			
			(Myrtan)		
			Tannic Acid		

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Differentiation of Some Natural Tannin Extracts by Fluorescence Analysis in Ultra - violet Light

HIROSHI MAEDA
(Central Customs Laboratory, 3-2 Kasumi - gaseki, Chiyoda - Ku, Tokyo)

(Received, Feb. 15, 1966)